



SHOWERS & EYEBATHS

Your Safety Demands Our Quality

Operating and Maintenance Instructions S300/E300 open eyebath.

The unit is manufactured according to ANSI Z358.1 - 2009. We recommend the unit be fed from potable water and where necessary a non-return check valve should be fitted to the supply pipe feeding the unit. The supply pipework should be sized at a minimum 1/2" diameter with a pressure of 2 bar g.

INSTALLATION INSTRUCTIONS

1. Position unit at correct height (approximately 1M from floor) in a suitable unobstructed well lit area.
2. Bolt to wall or bolt pedestal(if fitted) to floor
3. Connect water to inlet and flush out the eyebath to test for leaks and efficient operation.

MAINTENANCE INSTRUCTIONS

1. Ensure unit has a maintenance contract in place.
2. Test eyebath for good operation weekly as per ANSI.
3. Unit should be serviced and cleaned regularly at least every six months in order to ensure proper operation.
4. Clean the unit with non abrasive cleaners only. Replace any signs or tape when required.
5. Remove eyebath nozzles, inspect and disinfect every six months.
6. Replace eyebath nozzles annually
7. Check all fastenings every 6 months.
8. Report any non conformances.
9. Update service contract

OPERATING INSTRUCTIONS

1. Push plate to operate eyebath (or step on additional foot treadle if fitted).
2. To Turn off pull forward the push plate.
3. The valves remain open during the whole washing process.
4. Seek medical assistance after use.

Please note: If a scald valve is fitted, it is factory set at 38 deg C. It will release incoming water to the floor if the water temperature exceeds 38 deg C.

If a pressure reducing valve is fitted it is factory set at 2-3 bar to ensure correct flow from the unit.

Please do not adjust.

For any help or advise, please contact the maintenance department at Showers & Eyebaths Services Ltd.

World Class Manufacturers of Showers & Eyebaths

Showers & Eyebaths Services Ltd
Safety House, 23c Sandwash Close, Rainford Industrial Estate,
Rainford, Lancashire, WA11 8LY U.K.
Tel +44 (0) 1744 889677 Fax +44 (0) 1744 885663.
e-mail: sales@safety-showers.com www.safety-showers.com
Company Registration No. 3447603 Vat No: 703 9281 41



AMP 3 & 4 Framework Supplier